

M3 Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



EDAN

A world of potential

M3 Vital Signs Monitor

Proven to be a durable and efficient solution for outpatients and general wards, the M3 brings effective vital signs monitoring to most ambulatory and regular healthcares.

Measured Parameters: NIBP + SpO₂ + Infrared Ear Thermometer

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 14 hours of continuous work
- USB flash disk support
- HL7 support via XML files or MFM-CMS

30,000
NIBP Measurement Review

72 h
Trend Review

800
Alarm Review

- Barcode scanner support
- Bi-directional communication with MFM-CMS central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN/Wi-Fi support
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless Connectivity



Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp



Specification

| Configuration | SpO ₂ | NIBP | Infrared Ear Thermometer |
|---|------------------|------|--------------------------|
| SpO ₂ only | ✓ | — | — |
| NIBP only | — | ✓ | — |
| SpO ₂ +NIBP | ✓ | ✓ | — |
| NIBP+Infrared Ear Thermometer | — | ✓ | ✓ |
| SpO ₂ +NIBP+Infrared Ear Thermometer | ✓ | ✓ | ✓ |



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-M3-IET
V1.0-20200226