



SD3
Ultrasonic Pocket Doppler

SD3

Ultrasonic Pocket Doppler



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters:
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122 P.R. China
+86.755.26898326 | www.edan.com | info@edan.com.cn
U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com



ENG-OB & GYN-SD3-V2.5-20200101

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



SD3

Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by clinicians.



User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

Interchangeable Waterproof Probe


- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and larger patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- 4/5/8 MHz probe for blood flow detection
- Automatic power-off while putting the probe back to the rack (SD3 PLUS & SD3 PRO)

Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display

Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



Configurations

Model	SD3 VASCULAR	SD3 LITE	SD3	SD3 PLUS	SD3 PRO
FHR Display	×	×	√	√	√
Charger Stand	×	×	×	○	○
Record & Play	×	×	×	×	√
AA Batteries	○	○	○	×	×
Li-ion Battery	×	×	×	√	√

√ Standard × Not available ○ Optional



SD3 VASCULAR



SD3 LITE



SD3



SD3 PLUS



SD3 PRO