

# C3A/C6A Video Colposcope



C3A/C6A  
Video Colposcope



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Global Headquarters:**  
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,  
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122 P.R. China  
+86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com.cn](mailto:info@edan.com.cn)

**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)



© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

ENG-OB & GYN-C3A&C6A-V2.9-20200101



# C3A/C6A

Video Colposcope

With its bright LED cold lighting system, high-resolution camera and electronic green filter, C3A/C6A Video Colposcope realizes full-screen display of precise cervical examination, offering crucial information for clinical decision making to the gynecologist.

## Features

- LED cold lighting system reproduces an original view of tissue
- Advanced and fast auto-focusing system ensures clear image
- Remote capture control and user-friendly keyboard layout for easy operation
- Electronic green filter enhances vascular image display without compromise of illumination
- Acetic acid and iodine reaction test timer to capture the right and precise image
- Real-time display of magnification
- Powerful image display and data management software

## R-way®

The R-way® evaluation system reveals the variation of "color" on the cervix to interpret the performance of cervical diseases.

**CIN2+**

Saline application - 'R for Red'  
To observe the area with rich blood circulation, this reflects in red.

5% acetic acid-'W for White'  
To reveal increased nuclear activity, this reflects in white.

Green filter- 'A for Abnormal'  
To better observe abnormal vascular features, this defines as A.

Lugol's iodine solution- 'Y for Yellow'  
To observe fields of no iodine uptake, this reflects in yellow.

\* R-way® is not available for sale in the USA



## Patient management software

